



Hand-Held Shower

Model #: _____

Place of Purchase: _____

Serial #: _____

Phone #: _____

Date of Purchase: _____

Contact Person: _____

Setting up your equipment:

- If you use a shower chair, be sure the shower hose is long enough for use while seated in your chair.
- Follow manufacturer instructions to affix hand-held showerhead onto your existing shower arm.
- A diverter valve is included which controls where flow of water is directed – either through the hand-held showerhead or the faucet. **Tub or convention showerhead vs. hand-held showerhead in shower stalls.** You will choose the set-up at installation to suit your fixtures.

Using your equipment:

- Use control dial located on the showerhead to adjust water pressure.
- After use, secure hand-held shower onto the shower arm.
- Keep hose free of knots and kinks that would disrupt water flow.

Maintenance:

- Clean your hand-held shower with a disinfecting household cleaner as needed.
- The same cleaners used in the tub are appropriate.
- If a leak should occur and persist, call a plumber.

Safety Issues:

- Keep the hose free of knots, which could shorten the length of the hose and cause the user to trip or become entangled.
- If the shower is shared with those who stand to wash, courtesy suggests moving the diverter valve to the faucet position to prevent unexpected sprays from the hand-held showerhead.

Please contact your equipment provider for any questions, adjustments, or repairs.

Your local provider:



©Copyright The VGM Group 2002
B.L.U.E. Sheet photos, drawings, and text courtesy of
Golden Technologies, Invacare, Leisure-Lift, Permobil,
Puritan Bennett, Pride, Respironics, Ross,
Sammons Preston, Smith-Nephew, Sunrise Medical