



Tub Grab Bars

Model #: _____
Serial #: _____
Date of Purchase: _____

Place of Purchase: _____
Phone #: _____
Contact Person: _____

Setting up your equipment:

- Your provider may or may not install your tub grab bar for you. If they do not, they can refer you to a qualified professional to install it for you.
- Ask your provider to help you determine the proper positioning for the grab bar in relation to how you enter and exit the tub.

Using your equipment:

- When grab bar is securely fastened to the side of the tub, test it by applying pressure to all sides of the bar. If it is loose and/or appears to move slightly, tighten the grips to each side of the tub by turning the hand knob (Diagram 1). If you feel it is still not secure, contact your provider.


Maintenance & Replacement Items:

- Rubber accessories and grips may need to be replaced depending on wear.

Safety Issues:

- Test tub grab bar before initial use to be sure it is fastened securely to the tub side.
- Watch for cracks in your grab bar, and notify your provider of any significant and/or premature damage.

Diagram 1



Please contact your equipment provider for any questions, adjustments, or repairs.

Your local provider:



©Copyright The VGM Group 2002
B.L.U.E. Sheet photos, drawings, and text courtesy of
Golden Technologies, Invacare, Leisure-Lift, Permobil,
Puritan Bennett, Pride, Respironics, Ross,
Sammons Preston, Smith-Nephew, Sunrise Medical